



NOBA[®] DIALYSEUNTERLAGE

REF 003043

Product Description, Intended use, Application

The disposable product NOBA[®] DIALYSEUNTERLAGE is a 30 x 43 cm underlay consisting of a liquid-impermeable blue film back and a hydrophilic non-woven front with an inner cellulose layer. The underlay is closed all around the edges of the film to protect against liquid leakage. The cellulose layer absorbs liquids. The foil prevents the leakage of liquids into the environment. The dialysis pad supports the nursing staff in the clinic and in the ambulatory care of patients at home by collecting body fluids and soiling. The risk of contamination is thus reduced, and the spread of germs is counteracted.

Composition

Cellulose, polyethylene, polypropylene

Contra-Indications

The product should not be used in the case of a known allergy against the material.

Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Normative and Regulative Requirements, Common Standards

Medical device according to MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

Packaging

Primary packaging: PE-bag

Secondary packaging: carton made of cellulose

Symbols used in labelling

Explanations at www.nobamed.com



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

Storage and Transport

To be stored in a dry and dust-free environment

Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing, for reuse can decrease the product's performance significantly.

Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.