

Katheter-Set

REF 190100

Contents of this procedure pack:

1 x kidney-shaped tray, NOBAMED Paul Danz AG, MDC€

2 x gauze swab 7,5 x 7,5 cm, 8-ply, NOBAMED Paul Danz AG, MD (€

5 x gauze ball 3,0 cm, NOBAMED Paul Danz AG, MD < €

1 x NOBAWRAP[®] 60 x 60 cm wrapping drape, NOBAMED Paul Danz AG, MD C €

2 x NOBAGLOVE® LATEX size M, NOBAMED Paul Danz AG, MD C€

1 x forceps, blue, NOBAMED Paul Danz AG. MD(€

1 x NOBAWRAP® 60 x 60 cm, NOBAMED Paul Danz AG, MD C€

1 x plastic tray 125, NOBAMED Paul Danz AG, MDC €

This product is for single and temporary use on one patient and must be applied by specially trained medical personnel. Any information for use that is available for individual products must be observed!

Normative and Legal RequirementsProcedure Pack according to Art. 12 of the Medical Device Directive 93/42 / EEC.

The compatibility of the customized medical devices or other products has been assessed with due regard for the indications of

the respective manufacturers, and operations have been performed in accordance with the aforementioned indications.

The manufacturers' intended use is not changed.

The systems and procedure packs are packed in compliance with a validated procedure.

All activity is properly monitored and controlled as part of our QM system, in accordance with DIN EN ISO 13485:2016.

A validated sterilization process is used for sterile procedure packs, taking into account the manufacturers' instructions, so that product functionality and product safety of the components are maintained after a system or procedure pack has been compiled and sterilized.

Sterilization and packaging are subject to supervision of our Notified Body TÜV Süd PS 0123 pursuant to Annex V.

The product does not contain dangerous toxic substances according to REACH. It has CE marking and DIN EN ISO 15223-1 labels on all its packaging

Storage

To be stored in a dry and dust-free environment, protect from sunlight.

Symbols used in labellingExplanation at www.nobamed.com









