

### **NOBAMED Paul Danz AG**

## **NOBARIGID®-open**

#### **REF 943460**

## **Product Description, Intended use, Application**

NOBARIGID®-open (size S) is a ready-made frog-shaped aluminium splint with an inside padding on one side. It is used for immobilization of the finger. The product has to be bent into the desired finger position. To stabilize the splint, the protruding parts are bent to some degree around the finger. Then the splint is fixed by means of a bandage.

### Composition

Aluminium, polyester-based polyurethan foam

#### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

# Normative and Regulative Requirements, Common Standards

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

Primary packaging: PE bag

Secondary packaging: carton made of

cellulose

### Symbols used in labelling

Explanation at www.nobamed.com







### Storage and transportation

To be stored in a dry and dust-free environment

#### Single use device

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.