



## SLIPPER

**REF 944042**

### **Product Description, Intended use, Application**

Disposable shoe, size 42, in slip-on form for patients. It is packed in pairs and consists of blue non-woven and a sole made of EVAC. The product is worn in medical facilities such as hospitals and day clinics and prevents the spread of contamination.

### **Composition**

Polypropylene, EVAC (ethylene-vinyl acetate-copolymer)

### **Incident reporting**

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **Normative and Regulative Requirements, Common Standards**

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745. The product does not contain dangerous toxic substances according to REACH.

### **Packaging**

|            |                          |
|------------|--------------------------|
| Primary:   | PE pouch                 |
| Secondary: | carton made of cellulose |

### **Symbols used in labelling**

Explanations at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and DIN EN ISO 20417.

### **Storage and transport**

To be stored in a dry and dust-free environment

### **Single use device**

Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### **Disposal**

According to locally applicable legal regulations and standards of infection prophylaxis.