



## SUSPENSORIUM

REF 950006

### Product Description, Intended use, Application

The suspensory, size 6= S, is provided with a detachable bag. It is washable and is used for the fixation of the testicles after operations or in the case of oedema.

### Composition

Cotton, clip made of iron, buttons made of polypropylene

### Contraindications

The product should not be used in the case of a known allergy against the material.

### Incident reporting

According to MDR (EU) 2017/745, if serious incidents occur in relation to the device, they must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### Normative and Legal Requirements

Medical Device according to MDD 93/42/EEC, MDR (EU) 2017/745.

The product does not contain dangerous toxic substances according to REACH.

### Sizes

Size 4-6	child
Size 8	adolescent and adult
Size 10-12	seniors

### Packaging

Primary: folding box made of cellulose

Secondary: carton made of cellulose

### Symbols used in labelling

Explanation at [www.nobamed.com](http://www.nobamed.com)



Marking on all packaging levels with CE and according to DIN EN ISO 15223-1 and ISO 20417.

### Storage and Transport

Dry and dustfree

### Single use device

Intended for use on one patient. Reusing a single use medical device can lead to microbiological danger. Reprocessing for reuse can decrease the product's performance significantly.

### Disposal

According to locally applicable legal regulations and standards of infection prophylaxis.